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DISINFECTING AND CLEANING YOUR TONOMETER

(1) CLEANING

At the end of the day's use, detach the 5.5 gram weight and withdraw the plunger from the cylinder. Clean the plunger by wiping with a gauze sponge moistened with 70 to 95 percent alcohol. Thoroughly irrigate the plunger cylinder cavity with warm distilled water to dissolve any salt crystals which may have formed from tear fluid. Irrigate the cavity with 70 to 95 percent alcohol, and then allow tonometer to dry thoroughly. It is strongly recommended that the complete tonometer be stored in its closed case or in a sterilizing stand. Frequent inspection of the footplate and plunger is urged. Any nick or scratch caused by accidental mishandling must be corrected forthwith.

(2) DISINFECTION

Heat sterilization is recommended, using the Electric Tonometer Sterilizer which is designed to sterilize the plunger and footplate area with thermostatically controlled heat.

Cold sterilizing solutions which do not react with metals may be used in dilutions which will not be injurious or irritating to the patient's eyes. Consult manufacturers' recommendations before use.

Since the sterilizing solution may leave crystalline deposits or induce chemical attack on the stainless steel parts of the tonometer, it is important to clean the instrument after each such sterilization.

Any residue left on the tonometer may result in malfunction.